

ELISA Blocking Buffer



Catalog No.	Size
E104	1 packet

Description: Tris Buffered Saline with 1% BSA, pH 8.0. Sigma Chemical #T-6789, 1 packet

Application: Prepared Tris buffered saline with 1% BSA, pH 8.0, is suitable for use as a postcoat solution in ELISA.

Prepared Tris buffered saline with 1% BSA, pH 8.0, plus 0.05% Tween 20 (0.5 ml of 10% Tween 20 per 100 ml of postcoat) is suitable for use as sample and conjugate diluent.

Preparation: Dissolve contents of one packet in one liter deionized water to yield a 0.05 M Tris, 0.138 M NaCl, 0.0027 M KCl, pH 8.0 at 25 °C, with 1% BSA.

Storage: Packet can be stored at 2 – 8 °C before preparation.

Store prepared solution at 2 - 8 °C for up to one year. Discard if solution becomes cloudy.

For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Related Products

ELISA Starter Accessory Kit	E101
ELISA Starter Accessory Kit II	E103
TMB One Component Substrate	E102
Microtiter plate	E105
ELISA Wash Solution	E106
ELISA Coating Buffer	E107
10% Tween 20	E108
ELISA Stop Solution	E115