## Human IgG Heavy and Light Chain Cross-**Adsorbed Antibody**



Rabbit Polyclonal Conjugate DyLight® 488

Antigen Affinity Purified Catalog No. A80-218D2 A80-218D2-12 Lot No.

**APPLICATIONS** IHC, ICC, F, IF

SPECIES REACTIVITY Human. Minimum reactivity to mouse and rat

**AMOUNT** 1 ml

CONCENTRATION 0.5 mg/ml

STORAGE/SHELF LIFE 2 - 8°C / 1 year from date of receipt

PHYSICAL STATE Liquid

**BUFFER** Phosphate Buffered Saline (PBS) containing 0.2% BSA and 0.09% Sodium Azide

FLUOROPHORE/PROTEIN 7.3 **ISOTYPE** IgG **ORIGIN** USA

**PRODUCTION PROCEDURES** 

Antiserum was cross adsorbed using mouse and rat immunosorbents to remove cross reactive antibodies. The antibody to human IgG was isolated by affinity chromatography

using antigen coupled to agarose beads and conjugated to DyLight<sup>®</sup> 488.

Antibody concentration was determined by extinction coefficient: absorbance at 280 nm of

1.4 equals 1.0 mg of IgG.

By immunoelectrophoresis and ELISA this antibody reacts specifically with human IgG and with light chains common to other human immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. Less than 1% cross reactivity to mouse and

rat IgG was detected. This antibody may cross react with IgG from other species.

**APPLICATIONS** Centrifuge tube to remove product from lid. Optimal working dilutions should be determined

experimentally by the investigator. Prepare working dilution immediately before use.

Immunohistochemistry 1:50 - 1:500 Immunocytochemistry 1:50 - 1:500Flow Cytometry 1:50 - 1:200

**Immunofluorescence** 1:50 - 1:500

APPLICATION NOTES Not all listed applications have been specifically tested by our laboratory.

DyLight® 488 is excited at 493 (in PBS) and emits at 518 (in PBS).

DyLight<sup>®</sup> is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.

**ADDITIONAL INFO** https://www.fortislife.com/p/A80-218D2

Use the link above to view SDS, a current list of citations, and other product specific information.

This document certifies that this product has met all of the quality control standards defined by Bethyl Laboratories, Inc. Michael Spencer, PhD Date: July 19, 2022

> Phone: 800.338.9579 • Fax: 866.597.6105 • Web: www.fortislife.com Orders: orders@fortislife.com · Support: technical@fortislife.com