## Pig IgG-Fab Antibody

Goat Polyclonal Antiserum Catalog No. A100-113





Lot No.	A100-	113-1	LABORATORIES, INC
APPLICATIONS		IEP, DD	
SPECIES REACTIV	VITY	Pig	
ISOTYPE		IgG	
AMOUNT		2 ml	
STORAGE/SHELF	LIFE	2 - 8° C / 2 years from date of receipt	
PHYSICAL STATE	E	Liquid	
BUFFER		Serum containing 0.09% Sodium Azide	
ORIGIN		USA	
PRODUCTION PROCEDURES		By immunoelectrophoresis, the antiserum reacts specifically with Fab fragment of pig IgG and with light chains common to other pig immunoglobulins. No antibody was detected against non- immunoglobulin serum proteins.	
APPLICATIONS		Centrifuge tube to remove product from lid. Optimal working dilutions should be determined experimentally by the investigator. Prepare working dilution immediately before use.	
		IEP Neat	
		DD Neat	
APPLICATION NO	<b>LICATION NOTES</b> For use in precipitin gel reactions; such as, immunoelectrophoresis and double diffusion Ouchterlony.		and double diffusion
ADDITIONAL INF	FO	https://www.bethyl.com/product/A100-113 Use the link above to view SDS, a current list of citations, and other	product specific information.

This document certifies that this product has met all of the quality control standards defined by Bethyl Laboratories, Inc.Eric McIntush, PhD | Chief Scientific OfficerDate: December 3, 2018

For in vitro laboratory use only. Not for any clinical, therapeutic or diagnostic use in humans or animals. Not for human or animal consumption. This product may not be resold or modified for resale without the prior written approval of Bethyl Laboratories, Inc. The information provided in this data sheet is believed to be correct but does not purport to be all-inclusive and is intended to be used as a guide. Bethyl Laboratories, Inc. shall not be liable or responsible in any way for use of either this information or the material supplied. Disposal of hazardous material may be subject to federal, state or local laws or regulations.