## Monkey IgG Heavy and Light Chain Antibody



Goat Polyclonal Conjugate HRP

Antigen Affinity Purified Catalog No. A140-102P Lot No. A140-102P-31

**APPLICATIONS** WB, IHC, ICC, ELISA

SPECIES REACTIVITY Monkey **AMOUNT** 1 ml

CONCENTRATION 1 mg/ml

STORAGE/SHELF LIFE 2 - 8°C / 1 year from date of receipt

PHYSICAL STATE Liquid

Phosphate Buffered Saline (PBS) containing 0.2% BSA and 0.05% Pro-Clean 400 **BUFFER** 

**ISOTYPE** IgG USA **ORIGIN** 

The antibody was isolated by affinity chromatography using antigen coupled to agarose **PRODUCTION PROCEDURES** 

beads and conjugated to horseradish peroxidase (HRP).

Antibody concentration was determined by extinction coefficient prior to conjugation: absorbance at 280 nm of 1.4 equals 1.0 mg of IgG. Molar enzyme/antibody protein ratio is

4:1.

By immunoelectrophoresis and ELISA this antibody reacts specifically with monkey IgG and with light chains common to other monkey immunoglobulins. Species detected are Rhesus,

Baboon, Carb Green, Cyno, Chimp and African Green. Reaction to human was also detected. No antibody was detected against non-immunoglobulin serum proteins. This

antibody may cross react with IgG from other species.

**APPLICATIONS** Centrifuge tube to remove product from lid. Optimal working dilutions should be determined

experimentally by the investigator. Prepare working dilution immediately before use.

Western Blot 1:5,000 - 1:50,000

**Immunohistochemistry** 1:200 - 1:5,000

Immunocytochemistry 1:200 - 1:5,000

**ELISA** 1:10.000 - 1:100.000

**APPLICATION NOTES** Not all listed applications have been specifically tested by our laboratory.

**ADDITIONAL INFO** https://www.bethyl.com/product/A140-102P

Use the link above to view SDS, a current list of citations, and other product specific information.

This document certifies that this product has met all of the quality control standards defined by Bethyl Laboratories, Inc. Michael Spencer, PhD Date: June 16, 2022

> Phone: 800.338.9579 • Fax: 866.597.6105 • Web: www.bethyl.com Orders: orders@fortislife.com · Support: technical@fortislife.com