## Syrian Hamster IgG-heavy and light chain cross-adsorbed Antibody

Goat Polyclonal Conjugate Cy5

Antigen Affinity Purified
Catalog No. A140-201C5
Lot No. A140-201C5-3



**APPLICATIONS** IHC, ICC, IF, F

**REACTIVITY TESTED** Syrian Hamster. Minimum reactivity to bovine, chicken, horse, human, mouse, rabbit, rat & sheep.

**ISOTYPE** IgG

**AMOUNT** 1 ml at 0.5 mg/ml

STORAGE/SHELF LIFE 2 - 8° C / 1 year from date of receipt

PHYSICAL STATE Liquid

FLUOROPHORE/PROTEI 6.1

**BUFFER** Phosphate Buffered Saline (PBS) containing 0.2% BSA and 0.09% Sodium Azide

ORIGIN USA

PRODUCTION PROCEDURES

Antiserum was cross adsorbed using bovine, chicken, horse, human, mouse, rabbit, rat & sheep immunosorbents to remove cross reactive antibodies. The antibody to syrian hamster IgG was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to Cy5™.

Antibody concentration was determined by extinction coefficient: absorbance at 280 nm of 1.4 equals

1.0 mg of IgG.

By immunoelectrophoresis and ELISA this antibody reacts specifically with syrian hamster IgG and with light chains common to other syrian hamster immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. Less than 1% cross reactivity to bovine, chicken, horse, human, mouse, rabbit, rat & sheep IgG was detected. This antibody may cross react with IgG from other

species.

**APPLICATIONS** Centrifuge tube to remove product from lid. Optimal working dilutions should be determined

experimentally by the investigator. Prepare working dilution immediately before use.

Immunohistochemistry1:50 - 1:500Immunocytochemistry1:50 - 1:500Immunofluorescence1:50 - 1:500Flow Cytometry1:50 - 1:200

**APPLICATION NOTES** Not all listed applications have been specifically tested by our laboratory.

Cy5® is excited at 649 and emits at 670.

Cy® and CyDye® are registered trademarks of GE Healthcare.

ADDITIONAL INFO <a href="http://www.bethyl.com/product/A140-201C5">http://www.bethyl.com/product/A140-201C5</a>

Use the link above to view SDS, a current list of citations, and other product specific information.

