



Detection of human CLASP2 by western blot of immunoprecipitates. *Samples:* Whole cell lysate (1 mg for IP, 20% of IP loaded) from HeLa cells. *Antibodies:* Affinity purified rabbit anti-CLASP2 antibody A302-154A used for IP at 3 $\mu\text{g}/\text{mg}$ lysate. CLASP2 was also immunoprecipitated by rabbit anti-CLASP2 antibody A302-155A, which recognizes a downstream epitope. For blotting immunoprecipitated CLASP2, A302-155A was used at 1 $\mu\text{g}/\text{ml}$. *Detection:* Chemiluminescence with an exposure time 30 seconds.