Dog IgA Antibody

Rabbit Polyclonal

SPECIES REACTIVITY

Antiserum

Catalog No. A40–102 Lot No. A40–102–6

APPLICATIONS IEP, DD

ISOTYPE IgG

AMOUNT 2 ml

STORAGE/SHELF LIFE 2 - 8° C / 1 year from date of receipt

Dog

PHYSICAL STATE Liquid

BUFFER Serum containing 0.09% Sodium Azide

ORIGIN USA

PRODUCTION PROCEDURES

Antiserum was solid phase adsorbed to ensure class specificity.

By immunoelectrophoresis, the antiserum reacts specifically with dog IgA. No antibody was detected against immunoglobulin light chains or non-immunoglobulin serum proteins. This

antibody may cross react with IgA from other species.

APPLICATIONS Centrifuge tube to remove product from lid. Optimal working dilutions should be determined

experimentally by the investigator. Prepare working dilution immediately before use.

IEP Neat

DD Neat

APPLICATION NOTES For use in precipitin gel reactions; such as, immunoelectrophoresis and double diffusion

Ouchterlony.

ADDITIONAL INFO https://www.bethyl.com/product/A40-102

Use the link above to view SDS, a current list of citations, and other product specific information.

This document certifies that this product has met all of the quality control standards defined by Bethyl Laboratories, Inc.

Eric McIntush, PhD | Chief Scientific Officer

Date: December 3, 2018