Goat IgG-heavy and light chain Antibody

Rabbit Polyclonal

Antiserum

Catalog No. A50-100 Lot No. A50-100-14

APPLICATIONS IEP, DD SPECIES REACTIVITY Goat **ISOTYPE** IaG **AMOUNT** 2 ml

2 - 8° C / 2 years from date of receipt STORAGE/SHELF LIFE

PHYSICAL STATE Liquid

BUFFFR Serum containing 0.09% Sodium Azide

ORIGIN USA

PRODUCTION By immunoelectrophoresis, the antiserum reacts specifically with goat IgG and with light chains

PROCEDURES common to other goat immunoglobulins. No antibody was detected against non-

immunoglobulin serum proteins.

This antibody may cross react with IgG from other species.

APPLICATIONS Centrifuge tube to remove product from lid. Optimal working dilutions should be determined

experimentally by the investigator. Prepare working dilution immediately before use.

IEP Neat DD

APPLICATION NOTES For use in precipitin gel reactions; such as, immunoelectrophoresis and double diffusion

Neat

Ouchterlony.

https://www.bethyl.com/product/A50-100 **ADDITIONAL INFO**

Use the link above to view SDS, a current list of citations, and other product specific information.

This document certifies that this product has met all of the quality control standards defined by Bethyl Laboratories, Inc.

Eric McIntush, PhD | Chief Scientific Officer Date: December 3, 2018