Goat IgG Heavy and Light Chain Antibody



Donkey Polyclonal Conjugate DvLiaht® 488

Antigen Affinity Purified Catalog No. A50-101D2 A50-101D2-15 Lot No.

APPLICATIONS IHC, ICC, Flow Cyt, IF

SPECIES REACTIVITY Goat **AMOUNT** 1 ml

CONCENTRATION 0.5 mg/ml

STORAGE/SHELF LIFE 2 - 8°C / 1 year from date of receipt

PHYSICAL STATE Liquid

BUFFER Phosphate Buffered Saline (PBS) containing 0.09% Sodium Azide

FLUOROPHORE/PROTEIN 6.3 **ISOTYPE** IqG **ORIGIN** USA

PRODUCTION The antibody was isolated by affinity chromatography using antigen coupled to agarose **PROCEDURES**

beads and conjugated to DyLight® 488.

Antibody concentration was determined by extinction coefficient prior to conjugation:

absorbance at 280 nm of 1.4 equals 1.0 mg of IgG.

By immunoelectrophoresis and ELISA this antibody reacts specifically with goat IgG and with light chains common to other goat immunoglobulins. No antibody was detected against

non-immunoglobulin serum proteins.

This antibody may cross react with IgG from other species.

APPLICATIONS Centrifuge tube to remove product from lid. Optimal working dilutions should be determined

experimentally by the investigator. Prepare working dilution immediately before use.

Immunohistochemistry 1:50 - 1:500 **Immunocytochemistry** 1:50 - 1:500Flow Cytometry 1:50 - 1:200 Immunofluorescence 1:50 - 1:500

APPLICATION NOTES Not all listed applications have been specifically tested by our laboratory.

DyLight® 488 is excited at 493 (in PBS) and emits at 518 (in PBS).

DyLight® is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.

ADDITIONAL INFO https://www.fortislife.com/p/A50-101D2

Use the link above to view SDS, a current list of citations, and other product specific information.

This document certifies that this product has met all of the quality control standards defined by Bethyl Laboratories, Inc. Michael Spencer, PhD Date: December 6, 2022

> Phone: 800.338.9579 • Fax: 866.597.6105 • Web: www.fortislife.com Orders: orders@fortislife.com · Support: technical@fortislife.com