Horse IgG-heavy and light chain Antibody

Rabbit Polyclonal

Antiserum

AMOUNT

Catalog No. A70-118 Lot No. A70-118-3

APPLICATIONS IEP, DD
SPECIES REACTIVITY Horse
ISOTYPE IgG

STORAGE/SHELF LIFE 2 - 8° C / 2 years from date of receipt

2 ml

PHYSICAL STATE Liquid

BUFFER Serum containing 0.09% Sodium Azide

ORIGIN USA

PRODUCTION By immunoelectrophoresis, the antiserum reacts specifically with horse IgG and with light chains

PROCEDURES common to other horse immunoglobulins. No antibody was detected against non-

immunoglobulin serum proteins.

This antibody may cross react with IgG from other species.

APPLICATIONS Centrifuge tube to remove product from lid. Optimal working dilutions should be determined

experimentally by the investigator. Prepare working dilution immediately before use.

IEP Neat

DD Neat

APPLICATION NOTES For use in precipitin gel reactions; such as, immunoelectrophoresis and double diffusion

Ouchterlony.

ADDITIONAL INFO https://www.bethyl.com/product/A70-118

Use the link above to view SDS, a current list of citations, and other product specific information.

This document certifies that this product has met all of the quality control standards defined by Bethyl Laboratories, Inc.

Eric McIntush, PhD | Chief Scientific Officer



Date: December 3, 2018