Horse IgG-Fab Antibody

Rabbit Polyclonal

Antiserum

AMOUNT

Catalog No. A70-120 Lot No. A70-120-1

APPLICATIONS IEP, DD

SPECIES REACTIVITY Horse
ISOTYPE IgG

STORAGE/SHELF LIFE 2 – 8° C / 2 years from date of receipt

2 ml

PHYSICAL STATE Liquid

BUFFER Serum containing 0.09% Sodium Azide

ORIGIN USA

PRODUCTION PROCEDURES

Antiserum was solid phase adsorbed to ensure specificity.

By immunoelectrophoresis, the antibody reacts specifically with Fab fragment of horse IgG and with light chains common to other horse immunoglobulins. No antibody was detected against

non-immunoglobulin serum proteins.

This antibody may cross react with Fab Fragments from other species.

APPLICATIONS Centrifuge tube to remove product from lid. Optimal working dilutions should be determined

experimentally by the investigator. Prepare working dilution immediately before use.

IEP Neat

DD Neat

APPLICATION NOTES For use in precipitin gel reactions; such as, immunoelectrophoresis and double diffusion

Ouchterlony.

ADDITIONAL INFO https://www.bethyl.com/product/A70-120

Use the link above to view SDS, a current list of citations, and other product specific information.

This document certifies that this product has met all of the quality control standards defined by Bethyl Laboratories, Inc.

Eric McIntush, PhD | Chief Scientific Officer



Date: December 3, 2018