Human IaM Antibody

Rabbit Polyclonal		Conjugate	FITC		
Antigen Affinity Purified) · 5/ (
Catalog No. A80–101F Lot No. A80–101F–9				BETHYL	
	A80-10			LABORATORIES, INC	
APPLICATIONS		IHC, ICC, F, IF			
		Human			
AMOUNT					
CONCENTRATION					
STORAGE/SHELF LIFE		2 - 8° C / 1 year from date of receipt			
PHYSICAL STATE					
BUFFER		Phosphate Buffered Saline (PBS) containing 0.2% BSA and 0.09% Sodium Azide			
FLUOROPHORE/PROTEIN		3.6			
ISOTYPE		lgG			
ORIGIN		USA			
PRODUCTION PROCEDURES		Antiserum was solid phase adsorbed to ensure class specificity. The antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to fluorescein isothiocyanate (FITC).			
		Antibody concentration was determined by extinction coefficient prior to conjugation: absorbance at 280 nm of 1.4 equals 1.0 mg of IgG.			
		By immunoelectrophoresis and ELISA this antibody reacts specifically with human IgM. Cross reactivity with other immunoglobulins and light chains is less than 0.1%. This antibody may cross react with IgM from other species.			
APPLICATIONS		Centrifuge tube to remove product from lid. Optimal working dilutions should be determined experimentally by the investigator. Prepare working dilution immediately before use.			
		Immunohistochemistry 1:	50 - 1:500		
		Immunocytochemistry 1:	50 - 1:500		
		Flow Cytometry 1:	50 - 1:200		
		Immunofluorescence 1::	50 – 1:500		
APPLICATION NOTES		Not all listed applications ha	ve been specifically tested by our laboratory		
ADDITIONAL INFO https://www.bethyl.com/pr		https://www.bethyl.com/pro	oduct/A80–101F		
		Use the link above to view SI	DS, a current list of citations, and other prod	uct specific information.	

T

This document certifies that this product has met all of the quality control standards defined by Bethyl Laboratories, Inc. Eric McIntush, PhD | Chief Scientific Officer Date: August 20, 2020

Bethyl Laboratories, Inc. • 25043 West FM 1097 • Montgomery, TX 77356 • 800.338.9579 • 936.597.6111 • 866.597.6105 (FAX) • www.bethyl.com • technical@bethyl.com